

Listing Zeta 2017

Reference	Nom	Clones	Use	Vol.	Dilution
Z2064	Muscle Specific (HHF35	IVD	1ml	200
Z2066	Smooth Muscle (1A4	IVD	1ml	100
Z2068	AFP	Poly	IVD	1ml	100
Z2245	ALDH1A1	44	IVD	1ml	100
Z2289	ALK	1A4	IVD	1ml	30-50
Z2035	ALK/p80	SP8	IVD	1ml	100
Z2001L	AMACR	13H4	ASR	1ml	150
Z2005	AMACR+p63	13H4+4A4	ASR	0.5ml	50
Z2005L	AMACR+p63	13H4+4A4	ASR	1ml	50
Z2017	ACR+P63+HMW	3H4+4A4+34bE1	ASR	0.5ml	50
Z2017L	ACR+P63+HMW	3H4+4A4+34bE1	ASR	1ml	50
Z2139	Androgen Recepti	SP107	IVD	1ml	100
Z2246	Annexin A1	MRQ-3 (ZM51)	IVD	1ml	50-250
Z2220	Arginase-1	SP156	IVD	1ml	50
Z2283	ATRX	Poly	IVD	1ml	25-50
Z2120	β-Catenin	EP35	IVD	1ml	50-100
Z2318	BAP1	C4	IVD	1ml	100
Z2304	Cell Cocktail (p63+H	4A4 + 34BE12	IVD	1ml	50
Z2141	BCA-225	Cu-18	IVD	1ml	50
Z2022	BCL-2	124	IVD	1ml	100-200
Z2306	BCL-6	GI191E/A8	IVD	1ml	100-200
Z2260	BOB.1	SP92	IVD	1ml	50
Z2237	BRCA-1	MS110	IVD	1ml	50-100
Z2258	c-myc	EP121	IVD	1ml	50-100
Z2209	C3d	Poly	IVD	1ml	50-200
Z2210	C4d	Poly	IVD	1ml	50
Z2236	CA IX	EP161	IVD	1ml	50-100
Z2020	CA125	EP48	IVD	1ml	40-100
Z2072	CA19-9	121SLE	IVD	1ml	100-200
Z2286	Cadherin-17	SP183	IVD	1ml	50
Z2024	Calcitonin	SP17	IVD	1ml	50
Z2142	Caldesmon	E89	IVD	1ml	25-100
Z2122	Calponin-1	EP798Y	IVD	1ml	50-100
Z2026	Calretinin	SP13	IVD	1ml	50
Z2098	CD10	SP67	IVD	1ml	50
Z2030	CD117	YR145	IVD	1ml	50-250
Z2200	CD123	7G3	IVD	1ml	100
Z2251	CD13	SP187	IVD	1ml	50
Z2116	CD138	B-A38	IVD	1ml	20-100
Z2203	CD14	EP128	IVD	1ml	50-100
Z2028	CD15	MMA	IVD	1ml	100

Reference	Nom	Clones	Use	Vol.	Dilution
Z2259	CD16	SP175	IVD	1ml	50
Z2148	CD163	SP96	IVD	1ml	50
Z2185	CD19	SP110	IVD	1ml	100
Z2096	CD1a	O10	IVD	1ml	50-100
Z2308	CD2	EP222	IVD	1ml	25-100
Z2059	CD20	L26	IVD	1ml	100-300
Z2315	CD205	EP176	IVD	1ml	50-100
Z2123	CD21	EP64	IVD	1ml	50-200
Z2207	CD22	SP104	IVD	1ml	50
Z2053	CD23	SP23	IVD	1ml	50
Z2179	CD25	SP176	IVD	1ml	70
Z2033	CD3	Poly	IVD	1ml	50-200
Z2101	CD30	Ber-H2	IVD	1ml	50-100
Z2136	CD31	JC70A	IVD	1ml	50-100
Z2232	CD33	PWS44	IVD	1ml	20-80
Z2063	CD34	QBEnd/10	IVD	1ml	50-100
Z2241	CD35	EP197	IVD	1ml	50-100
Z2225	CD38	SP149	IVD	1ml	50
Z2273	CD4	EP204	IVD	1ml	50-100
Z2032	CD43	DF-T1	IVD	1ml	100
Z2080	CD44	SP37	IVD	1ml	50
Z2093	CD45/RB	2B11+PD7/26	IVD	1ml	100-200
Z2034	CD45RA	4KB5	IVD	1ml	50-100
Z2036	CD45RO	UCHL-1	IVD	1ml	50-100
Z2272	CD5	EP77	IVD	1ml	50-100
Z2081	CD54	23G12	RUO	1ml	10-20
Z2038	CD56	123C3.D5	IVD	1ml	100-200
Z2111	CD57	NK-1	IVD	1ml	100-200
Z2135	CD61	2f2	IVD	1ml	25-100
Z2071	CD68	KP1	IVD	1ml	100-200
Z2309	CD7	EP132	IVD	1ml	50-100
Z2268	CD71	MRQ-48 (ZM52)	IVD	1ml	50-100
Z2047	CD79a	SP18	IVD	1ml	100
Z2045	CD8	SP16	IVD	1ml	50
Z2250	CD99	EPR3097Y	IVD	1ml	50-250
Z2226	CDK4	DCS-31	IVD	1ml	50-100
Z2070	CDX-2	EPR2764Y	IVD	1ml	50-100
Z2100M	CEA	CEA31	IVD	1ml	100-200
Z2100P	CEA	Poly	IVD	1ml	100-200
Z2041	Chromogranin A	SP12	IVD	1ml	100
Z2102	CMV	8B1.2/1G5.2/2D4.2	ASR	1mg	100-200
Z2171	Collagen IV	PHM-12	IVD	1mg	50-100
Z2156	COX-2	SP21	IVD	1ml	50
Z2027	Cyclin D1	SP4	IVD	1ml	50
Z2307	Desmin	EP15	IVD	1ml	50-250
Z2279	Desmoglein-3	5G11	IVD	1ml	100
Z2013	DOG1	DOG1.1	IVD	1ml	100
Z2240	DPC4 (SMAD4)	B-8	IVD	1ml	50-100
Z2105	E-Cadherin	EP6	IVD	1ml	50-100

Reference	Nom	Clones	Use	Vol.	Dilution
Z2094	EBV-LMP	CS1-4	ASR	1ml	50-100
Z2037	EGFR	EP22	IVD	1ml	50-100
Z2288	EGFR (E746-A750 d	SP111	IVD	1ml	50
Z2287	EGFR (L858R)	SP125	IVD	1ml	50
Z2048	EMA	E29	IVD	1ml	100-200
Z2314	EP-CAM	Ber-EP4	IVD	1ml	50-200
Z2021	ER	ZR2	IVD	1ml	300
Z2021RS	ER	SP1	IVD	1ml	200
Z2217	ERCC1	SP68	IVD	1ml	50
Z2280	ERG	EP111	IVD	1ml	50-100
Z2162	ESA	MOC31	IVD	1ml	50-200
Z2323	EZH2	SP129	IVD	1ml	50
Z2107	Fa VIII-R Ag	Poly	IVD	1ml	10-50
Z2132	Factor XIIIa	EP3372	IVD	1ml	50-200
Z2079	Fasin	55K-2	IVD	1ml	50-250
Z2146	FLI-1	G146-22	IVD	1ml	25-50
Z2235	FoxP1	SP133	IVD	1ml	50
Z2301	FSH	EP257	IVD	1ml	50-100
Z2316	Galectin-3	B2C10	IVD	1ml	50
Z2303	Gastrin	Poly	IVD	1ml	50-250
Z2227	GATA3	L50-823	IVD	1ml	100
Z2073	GCDFP-15	EP1582Y	IVD	1ml	50-250
Z2255	GFAP	GA-5	IVD	1ml	200
Z2299	Glucogon	EP74	IVD	1ml	50-100
Z2264	GLUT1	EP141	IVD	1ml	50-100
Z2274	Glutamine Synthetas	GS-6	IVD	1ml	50-100
Z2103	Glypican-3	1G12	IVD	1ml	50-150
Z2145	Granzyme B	Poly	IVD	1ml	50
Z2300	Growth Hormone	EP267	IVD	1ml	50-100
Z2046	<i>H. Pylori</i>	Poly	ASR	1ml	100-200
Z2085	HBcAg	Poly	ASR	1ml	50-100
Z2331	HBsAg	A10F1	ASR	1ml	50-250
Z2108	HCG	CG04+CG05	IVD	1ml	50-100
Z2204	Hemoglobin A	EPR3608	IVD	1ml	50-250
Z2087	HepPar1	OCH1E5	IVD	1ml	100
Z2025	<i>Her-2/neu</i>	ZR5	IVD	1ml	200-300
Z2025E	<i>Her-2/neu</i>	EP3	IVD	1ml	100-200
Z2025RS	<i>Her-2/neu</i>	SP3	IVD	1ml	100
Z2201	HPV	K1H8	ASR	1ml	20
Z2215	HSV I	10A3	ASR	1ml	50-200
Z2216	HSV II	Poly	ASR	1ml	50-200
Z2211	IgA	Poly	IVD	1ml	50-200
Z2212	IgD	Poly	IVD	1ml	50-200
Z2213	IgG	Poly	IVD	1ml	50-200
Z2170	IgG4	MRQ-44 (ZM49)	IVD	1ml	100
Z2214	IgM	Poly	IVD	1ml	50-200
Z2284	IMP-3	EP286	IVD	1ml	50-100
Z2115	Inhibin- α	R1	IVD	1ml	10-100
Z2177	INI-1	25	IVD	1ml	50-100

Reference	Nom	Clones	Use	Vol.	Dilution
Z2298	Insulin	EP125	IVD	1ml	50-100
Z2297	Islet 1	EP283	IVD	1ml	50-100
Z2131	Kappa	L1C1	IVD	1ml	50-200
Z2042	Keratin 14	LL002	IVD	1ml	100-300
Z2254	Keratin 17	E3	IVD	1ml	100-200
Z2044	Keratin 18	DC10	IVD	1ml	200
Z2134	Keratin 19	BA17	IVD	1ml	100-200
Z2065	Keratin 20	Ks20.8	IVD	1ml	100-200
Z2257	Keratin 5	EP24	IVD	1ml	50-100
Z2133	Keratin 5/6	D5/16B4	IVD	1ml	100-200
Z2067	Keratin 7	OV-TL 12/30	IVD	1ml	50-100
Z2164	Keratin 8	35 β H11	IVD	1ml	40-100
Z2252	Keratin 8/18	B22.1 & B23.1	IVD	1ml	50-100
Z2018	Keratin CAM 5.2	CAM5.2	IVD	5ml	10
Z2018S	Keratin CAM 5.2	CAM5.2	IVD	1ml	10
Z2019	Keratin, HMW	34 β E12	IVD	1ml	100
Z2269	Keratin, HMW	AE3	IVD	1ml	100-200
Z2267	Keratin, LMW	AE1	IVD	1ml	100-200
Z2061	Keratin, Pan	AE1/AE3	IVD	1ml	100-200
Z2031	Ki-67	SP6	IVD	1ml	100
Z2305	Ki-67	MIB-1	IVD	1ml	150-300
Z2130	Lambda	Lamb14	IVD	1ml	100-500
Z2228	Langerin	12D6	IVD	1ml	25-100
Z2302	LH	Poly	IVD	1ml	50-250
Z2195	LMO2	SP51	IVD	1ml	50
Z2109	Lysozyme	Poly	IVD	1ml	50-200
Z2007	Mammaglobin	304-1A5	IVD	0.5ml	100
Z2007L	Mammaglobin	304-1A5	IVD	1ml	100
Z2009	Mammaglobin	31-A5	IVD	0.5ml	100
Z2009L	Mammaglobin	31-A5	IVD	1ml	100
Z2011	Mammaglobin Cock	304-1A5+31-A5	IVD	0.5ml	50
Z2011L	Mammaglobin Cock	304-1A5+31-A5	IVD	1ml	50
Z2295	MCM3	EP202	IVD	1ml	50-100
Z2189	MDM2	SPM14	IVD	1ml	50-100
Z2278	MDR-1	EP271	IVD	1ml	50-100
Z2052	Melan-A	A103	IVD	1ml	50-100
Z2088	Melanosome	HMB-45	IVD	1ml	100
Z2233	Mesothelial cell	HBME-1	IVD	1ml	25-100
Z2326	Mesothelin	EP140	IVD	1ml	50-100
Z2149	MGMT	MT3.1	IVD	1ml	20-40
Z2161	MiTF	C5/D5	IVD	1ml	50-200
Z2076	MLH-1	G168-728	IVD	1ml	100
Z2089	MMP-2	CA-4001	IVD	1ml	50
Z2090	MMP-9	EP127	RUO	1ml	50-100
Z2129	MSH-2	G219-1129	IVD	1ml	50-100
Z2281	MSH-6	EP49	IVD	1ml	50-100
Z2150	MUC-1	MRQ-17 (ZM35)	IVD	1ml	50-100
Z2151	MUC-2	Ccp58	IVD	1ml	100
Z2231	MUC-4	8G7	IVD	1ml	100

Reference	Nom	Clones	Use	Vol.	Dilution
Z2095	MUC-5AC	MRQ-19 (ZM37)	IVD	1ml	50-200
Z2163	MUC-6	MRQ-20 (ZM38)	IVD	1ml	50-100
Z2253	MUM1	MRQ-43 (ZR46)	IVD	1ml	50-200
Z2190	Myelin Basic Protein	Poly	IVD	1ml	50-200
Z2296	Myeloperoxidase	EP151	IVD	1ml	50-100
Z2261	MyoD1	EP212	IVD	1ml	25-100
Z2054	Myogenin	F5D	IVD	1ml	50
Z2126	Myoglobin	Poly	IVD	1ml	150-300
Z2229	Myosin, Smooth Muscle	SMMS-1	IVD	1ml	50-200
Z2294	Napsin A	ZM11	IVD	1ml	100-200
Z2178	NeuN	A60	IVD	1ml	50-100
Z2091	Neurofilament	2F11	IVD	1ml	50-200
Z2155	NGFR	MRQ-21 (ZM40)	IVD	1ml	50-200
Z2322	NKX2.2	EP336	IVD	1ml	50-250
Z2285	NKX3.1	EP356	IVD	1ml	50-100
Z2317	NRAS (Q61R)	EP174	IVD	1ml	50
Z2056	NSE	5E2	IVD	1ml	50-200
Z2310	OCT-2	EP115	IVD	1ml	50-100
Z2104	OCT-3/4	C-10	IVD	1ml	50-100
Z2249	Olig2	211F1.1	IVD	1ml	50-100
Z2192	p120 Catenin	SP63	IVD	1ml	50
Z2117	P16INK4 α	G175-405	RUO	1ml	100
Z2157	p21WAF1	DCS-60.2	IVD	1ml	50-200
Z2158	p27kip1	SX53G8	IVD	1ml	50-200
Z2004	P40	ZR8	IVD	1ml	100-200
Z2029	p53	DO-7	IVD	1ml	100-200
Z2173	p57	KP10	IVD	1ml	100-500
Z2003L	p63	4A4	RUO	1ml	50-100
Z2328	Parvalbumin	EP300	IVD	1ml	50-100
Z2292	PAX-2	EP235	IVD	1ml	20
Z2275	PAX-5	EP156	IVD	1ml	50-100
Z2312	PAX-6	EP341	IVD	1ml	50-100
Z2202	PAX-8	ZR1	IVD	1ml	50-100
Z2058	PCNA	PC10	IVD	1ml	100-200
Z2327	PD-1 (PDCD1)	EP239	IVD	1ml	50-100
Z2002	PD-L1	ZR3	IVD	1ml	200
Z2002S	PD-L1	ZR3	IVD	0.5ml	200
Z2118	Perforin	MRQ-23 (ZM44)	IVD	1ml	10-50
Z2193	PGP9.5	Poly	IVD	1ml	50-200
Z2208	PHH3	Poly	IVD	1ml	50-200
Z2313	PHOX2B	EP312	IVD	1ml	50-100
Z2218	PI3K P85	SP62	IVD	1ml	50
Z2043	PLAP	SP15	IVD	1ml	50
Z2270	PMS-2	EP51	IVD	1ml	20
Z2097	Podocalpin	D2-40	IVD	1ml	30-100
Z2023	PR	ZR4	IVD	1ml	400
Z2023RS	PR	SP2	IVD	1ml	400
Z2325	Prolactin	EP193	IVD	1ml	50-250
Z2006	Prostein (P501S)	ZR9	ASR	1ml	200

Reference	Nom	Clones	Use	Vol.	Dilution
Z2329	PSA	EP109	IVD	1ml	50-100
Z2140	PSAP	PASE/4LJ	IVD	1ml	100-200
Z2186	PSMA	SP29	IVD	1ml	50
Z2077	PTEN	6H2.1	IVD	1ml	50-100
Z2256	RCC	66.4.C2	IVD	1ml	100-200
Z2055	S-100	4C4.9	IVD	1ml	150-300
Z2277	S100P	EP234	IVD	1ml	50-100
Z2332	SALL4	EP299	IVD	1ml	50-100
Z2321	SATB2	EP281	IVD	1ml	50-100
Z2219	Smoothelin	R4A	IVD	1ml	50-100
Z2291	Insulin Receptor	EP149	IVD	1ml	50-100
Z2293	SOX-10	ZM10	IVD	1ml	100-200
Z2243	SOX-11	MRQ-58 (ZM50)	IVD	1ml	25-100
Z2198	SOX-2	SP76	IVD	1ml	50
Z2276	STAT6	EP325	IVD	1ml	50-100
Z2265	Stathmin	EP247	IVD	1ml	50-100
Z2196	Surfactant	1B9	IVD	1ml	50-100
Z2159	Survivin	SP79	IVD	1ml	50
Z2039	Synaptophysin	SP11	IVD	1ml	100
Z2324	TBX21 (T-bet)	EP263	IVD	1ml	50-100
Z2266	TCL-1	EP105	IVD	1ml	50-100
Z2271	TdT	EP266	IVD	1ml	50-100
Z2311	TFE3	EP285	IVD	1ml	50-100
Z2062	Thyroglobulin	2H11/6E1	IVD	1ml	50-100
Z2183	TIA-1	2G9A10F5	IVD	1ml	100
Z2248	TLE1	1F5	IVD	1ml	10-50
Z2175	Topoisomerase II α	Ki-S1	IVD	1ml	100-200
Z2166	TRAcP	9C5	IVD	1ml	50-200
Z2125	Tryptase	G3	IVD	1ml	50-200
Z2199	TS	TS106	IVD	1ml	50
Z2069	TTF-1	8G7G3/1	IVD	1ml	200
Z2074	Tyrosinase	T311	IVD	1ml	50-100
Z2083	Uroplakin III	SP73	IVD	1ml	50
Z2113	Villin	CWWB1	IVD	1ml	50-200
Z2051	Vimentin	SP20	IVD	1ml	100
Z2124	WT1	6F-H2	IVD	1ml	40-150
Z2114	ZAP-70	2F3.2	IVD	1ml	100-500